CRYSTAL YAG
HL 1364

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- Osteoartrosis
- Condrotipatías degenerativas

Therapeutic indications

The protocols developed, enable a quick resolution of the most common diseases in the fields of: orthopaedics, physiatry, dermatology, rheumatology and sports medicine.

The list of the main diseases which can be treated with it include:

- arthritis and degenerative of cartilage
- bursitis, synovitis, capsulitis, epicondylitis, shoulder impingement syndrome
- tendinitis and tenosynovitis
- oedema and trauma haematoma
- post-traumatic or overload pathology
- abductor syndromes, sprains, enthesitis, chondropathies kneecaps and sports traumatology.

More specifically, it is possible to treat:

1. Acute pathologies

- Tendinopathy (tendinitis, peritendinitis, tenosynovitis, tendinitis and calcific insertional, tendinosis, partial lesions)
- Trauma and post-traumatic oedema distorting
- Muscle injury
- Synovitis and bursitis
- Back pain and lumbago
- Fibromyalgia Syndrome
- Osteochondritis dissecans and osteochondral lesions

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Technical specifications	CRYSTAL YAG	CRYSTAL YAG DESK
Backlit LCD display	Graphic colour 8"	Graphic colour 8"
Peak power in continuous mode	25W	12W
Pulse adjustable	10 - 100%	10 - 100%
Setting emission	10 - 100%	10 - 100%
Emission	continuous and pulsed	continuous and pulsed
Emission frequency	1 - 1000 Hz	1 - 1000 Hz
Adjustable treatment time	Up to 99 minutes	Up to 99 minutes
Output channels	1 interchangeable	1 fixed
Stored protocols	33 man + 33 woman	33 man + 33 woman
Storable protocols on the internal memory	200	200
Storable protocols on the USB	•	•
Size	Trolley based 40 x 104 x 28 cm	39 x 17,5 x 28 cm
Weight	32 Kg	4,65 Kg

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CRYSTAL YAG and CRYSTAL YAG DESK are certified devices
KIWA CERMET ITALIA Spa CE 0476, according to Directive 93/42/EEC

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


The Physical Therapy by EME Physio

The laser power is featuring more and more the international medical scene, drawing a growing interest especially for its excellent effects that can be obtained according with its short time of application, thus resulting in a therapy that meets the high standards of modern life. Crystal YAG is the answer to the new needs of doctors and physiotherapists.

The therapeutic window

Crystal YAG supplies a wavelength of 1064 nm, which belongs to the "therapeutic window" of the lasers between 400 and 1600 nm. The light beam emitted by Crystal YAG, reaches up to 25W adjustable power, allowing a tissue penetration in depth. With these effects it is possible to make the photochemical bio-stimulation on the cell membrane and into the mitochondria acting deeply in the tissues and causing significant effects of metabolic stimulation, analgesia, anti-oedema, anti-inflammatory.

Crystal YAG makes it possible to safely treat all muscular-skeletal inflammation, superficial and deep, leading from the very first application to an intense pain reduction and a recovery in mobility.

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
Continuous mode and Duty Cycle

Crystal Yag can work in **continuous mode**, emitting a power up to 25W, or in **pulsed mode** Duty Cycle, allowing to switch from a time of laser emission to a non emission time, in order to reduce the "power effect" of the laser (i.e. photo-thermal effect which an increase in local temperature) and then allowing the cooling of the area treated between the different pulses. This makes it more effective and prominent the effect of bio-stimulation of the tissue.

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Effects on tissue

A significant stimulation of mitochondrial with an increase in:

- Production of ATP (photochemical effect: It accelerates the transformation of ADP into ATP, the electrolyte replacement between intra and extra cellular environments.)
- Activation of micro-circulation (photo-thermal effect: It activates the parts by changing the hydrostatic intracapillary pressure)
- Activating the lymphatic peristalsis (photo-mechanical effect: The lymphatic system is stimulated causing a rapid resorption of inflammatory agents with an effective action on fluid accumulation, too).
- Hyper-polarization of the nerve fibers membranes resulting in an analgesic effect (which raises the perceptual threshold of algotrope endings. It induces a significant increase in the release of endorphins.)

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The results

The Treatment with Crystal YAG provides:

- **Reduction of pain**, for some cases starting from the first treatment
- **Deep tissue action**
- **High bio-stimulating effect**
- **A winning combination**, which makes patients fully satisfied with the treatment.