



Welch Allyn Stethoscopes

Hear the Difference.



WelchAllyn[®]

Advancing Frontline Care[™]



Welch Allyn Stethoscopes

Hear the difference from the leader in diagnostic care

For almost 100 years, we have been focused on providing exceptional performance and value from every Welch Allyn device. Our complete lineup of stethoscopes are no exception.

From our innovative sensor-based Elite Electronic stethoscope, the Harvey DLX with trumpet brass construction for improved sound, to our well-known disposable line, every Welch Allyn product has been designed to provide you—the caregiver—with outstanding value and performance.

Electronic Stethoscopes

Hear more of what you need and less of what you don't

Elite Electronic with SensorTech™ Technology

Welch Allyn's sensor-based electronic stethoscope gives you precise sound pickup—minimizing vibration, resonance and sound loss. This makes it easier to detect low diastolic murmurs and high-pitched pulmonary sounds you might not hear with other stethoscopes.

Unlike traditional amplified stethoscopes that use a microphone in the head to collect the sound, the Elite Electronic places the acoustical sensor directly in contact with the diaphragm, allowing it to pick up even the faintest of sounds. Sound waves are then converted directly from the sensor to an electronic signal, minimizing unwanted ambient noise.

- 2-position filter switch lets you listen to high or low frequencies
- Adjustable binaurals and eartips ensure all-day comfort
- Long-life lithium ion battery lasts up to 140 hours
- Fits comfortably around your neck or in your pocket
- 1-year warranty

The Analyzer—Phonocardiogram and ECG right on your PC
Record, playback, analyze, store, and share files right from your PC.

The Distributor—Share what you hear

Share what you hear with up to five people at the same time.



ELITE ELECTRONIC
STETHOSCOPE 5079-405



ANALYZER SOFTWARE
5079-430



DISTRIBUTOR 5079-401

Cardiology Grade Stethoscopes

Developed during a 30-year association with Dr. W. Proctor Harvey, MD, a recognized "Master" by the American College of Cardiology.

Harvey DLX™

With trumpet brass construction, the Harvey DLX stands out among all other stethoscopes. Trumpet brass construction translates to hearing more of the important cardiac, pulmonary and vascular sounds without relying on difficult "variable pressure" techniques.

- Full-range cardiology/pulmonology stethoscope
- Available in double- or triple-head designs
- Flat or corrugated diaphragms
- Pediatric chestpiece available
- Chrome-plated brass chestpiece and binaurals
- Rotatable binaurals and interchangeable comfort sealing eartips
- Dual-lumen tubing creates individual, uninterrupted sound channels to each ear
- Latex free
- 10-year warranty

Harvey™ Elite®

The Harvey Elite offers updated styling with the superior comfort and auscultation capabilities cardiologists have come to rely on. Its stainless steel head transmits sounds crisply and clearly.

- Full-range cardiology/pulmonology stethoscope
- Double head (bell and diaphragm)
- Pediatric chestpiece available
- Stainless-steel chestpiece with non-chill rim
- Rotatable binaurals and interchangeable comfort sealing eartips
- Dual-lumen tubing creates individual, uninterrupted sound channels to each ear
- Flat or corrugated diaphragms
- Latex free
- 10-year warranty

NEW: 28" LENGTH AVAILABLE



HARVEY DLX DOUBLE-HEAD STETHOSCOPE 5079-325



HARVEY ELITE DOUBLE-HEAD STETHOSCOPE 5079-125

Professional Grade Stethoscopes

Excellent acoustic clarity and durability

Welch Allyn's Professional and Lightweight range of stethoscopes provide excellent sound quality, durability and comfort for the practicing professional. They are lightweight and can be custom fitted via adjustable binaurals and eartips—perfect for listening to cardiac, pulmonary and vascular sounds.

Professional Adult

- Excellent for general auscultation and blood pressure
- Double head design
- Flat diaphragm
- Adult and Pediatric versions available
- Stainless steel chestpiece and binaurals
- Rotatable binaurals and interchangeable comfort sealing eartips
- Latex free
- 2-year warranty

Lightweight

- Excellent for general auscultation and blood pressure
- Chrome-plated chestpiece for durability
- Rotatable binaurals
- Double head design
- Flat diaphragm
- Adult version only
- Latex free
- 1-year warranty



PROFESSIONAL ADULT STETHOSCOPE 5079-135



LIGHTWEIGHT STETHOSCOPE 5079-73

Disposable Stethoscopes (Single Patient Use)

The best protection against cross-contamination

Welch Allyn Uniscope is a single-patient-use stethoscope for use in situations where cross-contamination is a concern. It is lightweight and offers excellent acoustics. It is available as an Adult version (Yellow) and Pediatric version (Blue). The Uniscope allows accurate auscultation without compromising primary stethoscopes or patient safety.

Uniscope

- Single-patient-use (disposable) stethoscope
- Excellent for general auscultation
- Adult and Pediatric versions available
- Plastic chestpiece and binaurals
- Single-head design
- Flat diaphragm
- Latex free



UNISCOPE 17461

Welch Allyn Stethoscope Specifications

	ELITE ELECTRONIC STETHOSCOPE	CARDIOLOGY GRADE		PROFESSIONAL GRADE			DISPOSABLE
PRODUCT DESCRIPTION	ELITE ELECTRONIC	HARVEY DLX (28")	HARVEY ELITE DOUBLE HEAD (28") (22")	PROFESSIONAL ADULT	PROFESSIONAL PEDIATRIC	LIGHTWEIGHT	UNISCOPE
PRODUCT NUMBERS	5079-405 5079-430 (Analyzer) 5079-401 (Distributor)	5079-321 Triple, Black 5079-322 Triple, Burgundy 5079-323 Triple, Navy 5079-324 Triple, Green 5079-325 Double, Black 5079-326 Double, Burgundy 5079-327 Double, Navy 5079-328 Double, Green	5079-122 Elite, Black (22") 5079-125 Elite, Black 5079-270 Elite, Burgundy 5079-271 Elite, Navy 5079-284 Elite, Green 5079-125P Elite, Ped Blk	5079-135 Black 5079-137 Navy 5079-289 Blue 5079-139 Burgundy 5079-285 Forest	5079-145 Black 5079-147 Navy 5079-149 Burgundy 5079-287 Forest	5079-73 Velvet Black 5079-74 Poppy Red 5079-75 Dawn Gray 5079-76 Misty Blue	17461 Adult, Yellow 17462P Pediatric, Blue
PRODUCT TYPE	Electronic Stethoscope	Double/Triple Head	Double Head	Double Head	Double Head	Double Head	Single Head
HEAD DIAMETER INCHES (mm)	1.5 (38.5)	1.9 (48.3)	1.9 (48.3)	1.75 (45)	1.2 (30)	1.9 (48.3)	1.9 (49)
OVERALL LENGTH INCHES (mm)	32.5 (825.5)	25 (635)	25 (635)	28 (711)	28 (711)	32 (812.5)	29.5 (750) 27 (686)
TUBING LENGTH INCHES (mm)	21.5 (546.1)	17 (431.8)	17 (431.8)	21.5 (546.1)	21.5 (546.1)	24 (609.6)	19 (483)
LATEX FREE (Y/N)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
CHESTPIECE MATERIAL	Plastic	Brass	Stainless Steel	Stainless Steel	Stainless Steel	Chrome Plated Zinc	ABS
DIAPHRAGM MATERIAL	Plastic	Polycarbonate	Polycarbonate	Clear Polycarbonate	Iten Resiten G-10/FR4T Glass Epoxy Sheet	PVC	Polypropylene
DIAPHRAGM STYLE	Flat	Flat or Corrugated	Flat or Corrugated	Flat	Flat	Flat	Flat
Adjustable Binaural (Y/N)	Yes	Yes	Yes	Yes	Yes	Yes	No
EARTIPS (SOFT/HARD)	Soft	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Hard	Hard
TYPE OF TUBING (SINGLE/DUAL LUMEN)	Single Lumen	Dual Lumen	Dual Lumen	Single Lumen	Single Lumen	Single Lumen	Single Lumen
ANTI-CHILL RING (Y/N)	No	Yes (Brass)	Yes	Yes	Yes	Yes	No
SHIPS ASSEMBLED (Y/N)	Yes	Yes	Yes	Yes	Yes	No	Yes
PACKAGING	Printed Box	Printed Box	Printed Box	Printed Box	Printed Box	Plain White Box 1-piece	Plastic Bag
DESCRIPTION	Meets all diagnostic needs of the demanding professional. Provides the user with the ability to filter and amplify heart and lung sounds during auscultation. External 2.5 mm jack can be used to connect the stethoscope to an external device such as the Analyzer or the Distributor. Excellent tool for those experiencing hearing loss. Uses one 3 V lithium ion battery. Latex free.	Designed to meet the demanding needs of cardiologists, this stethoscope has fully adjustable binaurals and triple leaf spring for comfort. The head is made of trumpet brass because of its ability to transmit sounds accurately. The dual-lumen tubing is latex free.		Meets all diagnostic needs, from detecting Korotkoff sounds of blood pressure to analyzing subtle heart murmurs. Stainless-steel chestpiece intensifies body sounds transmitted through tissues. Comfort sealing eartips depress and conform to the shape of ear canals for soft yet airtight fit. Latex free.		Quality design for blood pressure measurement and general auscultation. Lightweight and flexible single-channel tubing provides comfortable use over extended periods. Adjustable chrome-plated binaurals. Latex-free.	The Uniscope is a single-patient-use stethoscope for situations where cross-contamination is a concern.

4341 State Street Road, PO Box 220, Skaneateles Falls, NY 13153-0220 USA
(p) 1.800.535.6663 (f) 315.685.3361 www.welchallyn.com

WelchAllyn[®]

Advancing Frontline Care™